

Smokerlyzer®

COmpact™ Smokerlyzer® ca. ppm & %COHb

babyCOmpact Smokerlyzer®

special
med®
Specialmed GmbH
Roßberg 2
D-84164 Lengenthal
Tel.: +49 8731 3264130
Fax: +49 8731 3264930
E-Mail: info@specialmed.de
Page: www.specialmed.de

	Micro+/piCO+ Erwachsene			Micro+/piCO+ Jugendliche/Schwangere			
	%COHb	COppm		COppm	%COHb	%FCOHb	
7	10.23	60		60	10.23	Mutter	Partner (Passiv Raucher)
	10.07	59		59	10.07		
	9.91	58		58	9.91		
	9.75	57		57	9.75		
	9.59	56		56	9.59		
	9.43	55		55	9.43		
	9.27	54		54	9.27		
	9.11	53		53	9.11		
	8.95	52		52	8.95		
	8.79	51		51	8.79		
6	8.63	50		50	8.63	5.66	2.64
	8.47	49		49	8.47		
	8.31	48		48	8.31		
	8.15	47		47	8.15		
	7.99	46		46	7.99		
	7.83	45		45	7.83		
	7.67	44		44	7.67		
	7.51	43		43	7.51		
	7.35	42		42	7.35		
	7.19	41		41	7.19		
5	7.03	40		40	7.03		
	6.87	39		39	6.87		
	6.71	38		38	6.71		
	6.55	37		37	6.55		
	6.39	36		36	6.39		
	6.23	35		35	6.23		
	6.07	34		34	6.07		
	5.91	33		33	5.91		
	5.75	32		32	5.75		
	5.59	31		31	5.59		
4	5.43	30		30	5.43		
	5.27	29		29	5.27		
	5.11	28		28	5.11		
	4.95	27		27	4.95		
	4.79	26		26	4.79		
	4.63	25		25	4.63		
	4.47	24		24	4.47		
	4.31	23		23	4.31		
	4.15	22		22	4.15		
	3.99	21		21	3.99		
3	3.83	20		20	3.83	3.40	1.49
	3.67	19		19	3.67		
	3.51	18		18	3.51		
	3.35	17		17	3.35		
	3.19	16		16	3.19		
	3.03	15		15	3.03		
	2.87	14		14	2.87		
	2.71	13		13	2.71		
	2.55	12		12	2.55		
	2.39	11		11	2.39		
2	2.23	10		10	2.23	2.38	1.20
	2.07	09		09	2.07		
	1.91	08		08	1.91		
	1.75	07		07	1.75		
	1.59	06		06	1.59		
	1.43	05		05	1.43		
	1.27	04		04	1.27		
	1.11	03		03	1.11		
	0.95	02		02	0.95		
	0.79	01		01	0.79		

www.bedfont.com

